

**BESTEC**

**HZB** Helmholtz  
Zentrum Berlin

Freie Universität Berlin



**MEDSI** PARIS  
MECHANICAL ENGINEERING DESIGN OF SYNCHROTRON  
RADIATION EQUIPMENT AND INSTRUMENTATION  
2018

Visit the BESTEC booth (no.16) for information about more spectrometers.

INSTRUMENTATION FOR SYNCHROTRON RADIATION

# High-efficiency and compact von Hamos spectrometer for the soft X-ray range

Filip Fuchs<sup>1</sup>, Tatjana Giebel<sup>1</sup>, Franz Schäfers<sup>2</sup>, Burkhard Langer<sup>3</sup>, Eckart Rüh<sup>3</sup>

<sup>1</sup>Bestec GmbH, Am Studio 2b, D-12489 Berlin, Germany,

<sup>2</sup>Helmholtz-Zentrum Berlin (HZB) BESSY II, Albert-Einstein-Str. 15, D-12489 Berlin, Germany,

<sup>3</sup>Freie Universität Berlin, Takustr. 3, D-14195 Berlin, Germany


## Introduction & Motivation

Samples illuminated by pulsed soft X-rays with a large number of photons per pulse (XFEL, Laser-plasma sources) often restrict analytical methods to ph-in-ph-out techniques, since ph-in-e-out techniques can suffer severely from space charge.

We introduce a compact and large angle acceptance soft X-ray spectrometer for an energy range of 1-3 keV suitable for experiments under the above conditions.

## Concept


The spectrometer uses a von Hamos optical design and covers a simultaneous energy window of >300 eV over the whole energy range.



Principle of the von Hamos design, where a thin crystal is bent to a cylindrical surface [1,2].

Performance defining criteria:


- Crystal lattice constant
- Bragg reflex line width
- Crystal radius
- Spatial detector resolution
- Crystal figure error
- Alignment & stability





→ Choice of crystal: Mica, 40 mm radius

→ Target effective tangent error: <100 μrad


Mica crystal structure.



Material	Lattice constant
Graphite	0.335 nm
Quartz	0.425 nm
Mica	0.99 nm
KAP	1.332 nm




## Layout




Using a high speed camera (500.000 frames/s at 1024x16 pixels) up to 400 spectra can be recorded within a macro bunch (10 Hz) of the European XFEL.

## Performance test results



Internal mechanics of the von Hamos spectrometer.



Diffraction from a thin mica crystal measured at the Optics beamline at BESSY/HZB.

## Specifications

Parameter	Typical value
Energy range	800 – 2300 eV
Energy resolution	0.6 – 2.5 eV
Energy spectrum width	~300 eV (simultaneous)
Azimuthal acceptance	>100°
Crystal	Mica
Crystal radius / length	40mm / 25 mm
Detector	MCP
Mechanical repeatability	< 0.5 μrad / < 0.1 μm
Mechanical stability	< 0.2 μrad / < 0.1 μm
Dimensions without camera	DN250 CF x 650 mm
Dimensions with camera	DN250 CF x 1800 mm

**BESTEC** GmbH  
Am Studio 2b, 12489 Berlin, Germany  
www.bestec.de  
phone: +49 30 677 4376  
fax: +49 30 677 5718



 Bundesministerium  
für Bildung  
und Forschung

## References

- [1] L. von Hamos, *Naturwissenschaften* **20**, 705 (1932).
- [2] A. P. Shevelko *et al.*, *Review of Scientific Instruments* **73**, 3458 (2002).

Financial support by the BMBF (contract no.: 05K13KE2) is gratefully acknowledged.