



Beam Loss Position Monitoring with Optical Fibres at DELTA

Frank Rüdiger

M. Körfer (DESY)
W. Göttmann (HMI Berlin)
G. Schmidt (DELTA)
K. Wille (DELTA)

24. June 2008




Table of Content

1. Introduction and Motivation
2. Optical Time Domain Reflectometry
3. Transmission Measurement with Optical Fibres
4. Detection of Cerenkov-Radiation
5. Summary

1. Overview of DELTA

- 1.5 GeV
- 130 mA
- lifetime: 10h
- 115 m



1. Radiation Sources at DELTA

- The DELTA vacuum chamber is made up of 3 mm V4A steel
- The synchrotron radiation is almost completely absorbed inside the chamber wall
- Beam loss electrons colliding with the chamber generate electro-magnetic cascades. Even at small incident angles a significant amount of the shower particles can leave the vacuum chamber



- **1.5 GeV beam loss electrons are the main source of ionising radiation at DELTA**

1. Challenges for Dosimetry Systems / Advantages of Fibre Optical Systems

- Measurements in narrow spaces
→ small fibre diameter
 - Surveillance of large distances
→ complete surveillance of
several kilometres
 - High-dose measurements
→ up to 1000 Gray
 - Evaluation during beam operation
→ possible within a few minutes
- System with two nanoseconds
time resolution is available




1. Radiation Induced Attenuation of Optical Fibres

- Chemical bonds are split up in the fibre by exposure to radiation (radiolysis)



- The generated defects are called “colour-centres”
- Transitions between the generated states
→ attenuation of injected light intensity

1. Radiation Induced Attenuation of Optical Fibres



- The sensitivity of a fibre optic radiation sensor system can be chosen by wavelength selection

1. Characteristics of a Fibre Optic Radiation Sensor System


- Linearity between dose α_D and attenuation D

$$\alpha_D = c(\lambda) \cdot D^f \quad [\text{dB/m}] \quad f \cong 1$$

valid for a dose range of 10^{-3} to 10^3 Gy

- Annealing: loss of accumulated dose information due to regeneration processes of the colour centres.
Enhancement of dose information lifetime by suitable fibre doping
- The attenuation is independent of dose-rate, temperature and light intensity

2. Optical Time Domain Reflectometry (OTDR)



- Used for radiation dose measurement around the complete length of the storage ring. Resolution about 3 Gray
- Measurement of the Rayleigh backscattered part of the injected light intensity
- The time dependent signal is converted to the position of beam loss; spatial resolution ca. 60 cm

2. OTDR: Hardware



- Tektronix TFP2A
- Selectable wavelength (850 & 1300nm)
- Pulse lengths 1ns, 3 ns, 8 ns, ...

Used Fibres:

- Multi-mode fibres with core diameter of 50 μm
- Germanium doped and co-doped with Phosphorus


3. Transmission Measurement: Motivation

Radiation effects of the permanent magnet undulator U55



- Magnets are composed of Neodymium-Iron-Boron alloy
- High radiation doses result in an irreversible demagnetization (based on investigations at the ESRF)
- Limiting dose value: 60 kGy
 - limiting dose rate value for the U55:
0.8 Gy/h
 - planned frequent injection mode will increase the radiation dose
 - **permanent dose surveillance of the U55 needed**



3. Transmission Measurement: Principle and Setup




3. Transmission Measurement: Hardware



- Powermeter: Newport Type: PTS-FOPM
- Detectable intensity range: 1pW – 2W
- Used wavelength: 660 nm
- The system has been calibrated using thermoluminescence dosimeter-rods
- accuracy is about 30 %, similar to the TTF-system at DESY

3. Results: Measured Fibre Intensity and Beam Current

Power Loop 1 @ 1.5 GeV Multibunch





3. Results: Dose Rates During Standard User Operation

- Measurement during 60 hours at 1.5 GeV during standard user operation; beam lifetime: 4h @ 100 mA
- Gap is opened to maximum during beam injection mode

	U55 entrance (upper side)	U55 entrance (lower side)
Dose rate [Gy/h]	0.81 ± 0.23	0.16 ± 0.05


- Limiting value already reached at the upper side
- The dose at the lower side is by a factor 5 smaller

3. Results: Dose Rates during Frequent Injection Mode

- Measurement during 2 hours at 1.5 GeV, frequent injection






	U55 entrance (vacuum chamber)
Dose rate [Gy/h]	4.23 ± 1.23

- the limiting value is exceeded by factor 5



- First measure: Optimisation of the lead shielding in front of the U55
→ dose rate lowered by factor 2 → sufficient for standard beam operation
→ **redesign of the lead shielding until start of frequent injection mode**

4. Cerenkov-Light Detector: runtime measurement to localise beam losses



Frank Rüdiger



4. Cerenkov-Light Detector: Principle of Measurement

- Runtime-measurement to localise beam losses
→ detector is connected to the upstream-side
- **Time resolution: 2 ns** given by sampling rate of the ADC (1 GS/s)
- **Spatial resolution: 0.24 m** in longitudinal direction
- Used wavelength range: 500 nm – 650 nm (maximum at 550 nm)
- Multi-mode-step-index fibres (core diameter: 300 μm) consisting of undoped silicon dioxide with high content of OH^- ions



5. Summary:

- Fibre optic radiation sensor systems are **well suited** for accelerators:
 - usable in narrow spaces
 - dose range is up to 1000 Gy
 - evaluation during beam operation
- **OTDR:**
 - used for dose surveillance of the complete Delta vacuum chamber
 - dose resolution: 3 Gy
- **Transmission measurement:**
 - used for dose surveillance of the U55 permanent magnet undulator
 - dose resolution: 60 mGy
- **Cerenkov-light detector:**
 - system has been installed and functionality has been proven
 - will be used for increase of injection efficiency
 - real-time beam loss position monitoring with single bunch resolution of 2 ns



Acknowledgement

- Dr. Markus Körfer (DESY)
- Walter Göttmann (HMI Berlin)
- Prof. Dr. Wille (DELTA)
- Prof. Dr. Weis (DELTA)
- Dr. Gerald Schmidt (DELTA)
- Holger Huck (DELTA)



References:

- Frank Rüdiger, diploma thesis, Aufbau und Einsatz von Glasfaser-Dosimetriesystemen an der Speicherringanlage DELTA, 2008
- Manuel Benna, diploma thesis, Strahlverlustmessung und Dosimetrie am DELTA, 2006
- S. Girard, A. Boukenter, et al., Properties of phosphorus-related defects induced by γ -rays and pulsed X-ray irradiation in germanosilicate optical fibres, 2003
- D. Griscom, E. Friebele, Fundamental radiation-induced defect centers in synthetic fused silicas: Atomic chlorine, delocalized E' centers and a triplet state, 1986
- H. Henschel, M. Körfer, J. Kuhnenn, U. Weiland, F. Wulf, Fibre optic radiation sensor systems for particle accelerator, Nuclear Instruments and Methods in Physics Research A 526, 2004
- H. Henschel, O. Köhn, H.U. Schmidt, Optical Fibres as radiation dosimeters, Nuclear Instruments and Methods in Physics Research B 69, 1992
- G. Schmidt et al., Optical Fibre Beam Loss Monitors for Storage Rings at DELTA, Proceedings of EPAC, 2002
- W.V. Hassenzahl, T.M. Jenkins et al., An Assessment of the Radiation on Permanent-Magnet Material in the ALS Insertion Devices, Nuclear Instruments and Methods in Physics Research A291, 1990
- G. Schmidt et al., Proposal for a Frequent Injection Mode at DELTA, Proceedings of EPAC, 2004
- H. Henschel et al, Optical Fibre Dosimeter For SASE FEL Undulators, DIPAC, 2003
- W. Göttmann, M. Körfer, J. Kuhnenn, Beam Loss Position Monitoring using Cerenkov Radiation on Optical Fibers, Proceedings of DIPAC, 2005