

Gas Sheet Ionization Diagnostic for Transverse Profile Measurement


N. Burger

RadiaBeam


NAPAC 2022

Albuquerque, NM


- Motivation
- Gas Sheet Ionization Monitor
- Gas Sheet Generation – Design
- Ionization Dynamics – Facet II
- Ion Microscope – Design
- Gas Sheet Generation – Validation
- Profile Reconstruction
- Commissioning Status
- Summary




- High-intensity beams present unique challenges in transverse profile diagnostics
- Single-shot, regenerative diagnostic
- “Gas sheet ionization monitor”
 - Generate a gas sheet, or “curtain” with nozzles and beam skimmers
 - 45° curtain for both axes projections in a single shot
 - Main particle beam ionizes neutral gas
 - Ions imaged by an “ion microscope”
 - Resolve high-intensity beams



Beam-gas interaction




- Gas sheet generation, vacuum considerations
 - Externally triggered, remote controlled Piezo-valve for tunable gas pulse
 - Modular system accepts many skimmer sizes
 - High-capacity pumping
- Ionization and extraction/imaging
 - Tunable electrostatic lens array
 - Microchannel plate / Phosphor screen
- Data analysis and ML beam reconstruction




Optimization of gas density profile: (Molflow+)

- Example to study effect of skimmer size/shape/locations
- Compact design to maintain distribution at IP



Differential pumping: (Molflow+)

- Strict UHV requirements determine pump speed and placement





Vacuum levels near IP (mbar)

	Min	Max	Av
Top	1,66E-09	4,49E-09	3,14E-09
Bottom	1,65E-09	4,13E-09	2,64E-09
Left	2,09E-09	4,54E-09	3,20E-09
Right	1,73E-09	4,50E-09	3,19E-09

Configuration 1: drive/witness

$Q = 0.5 \text{ nC}$
 $\sigma_x = 5 \mu\text{m}$,
 $\sigma_y = 7.5 \mu\text{m}$,
 $\sigma_z = 14 \mu\text{m}$
 Peak field
 $\sim 20 \text{ GV/m}$





N_2 gas sheet
 Density = 10^{19} m^{-3}
 Thickness = $150 \mu\text{m}$
 $Q_{\text{ions}} = 0.2 \text{ fC}$ (10^3 ions)

- Lower end of acceptable statistics
- High gain MCP
- Thicker gas sheet

Configuration 2: max compression

$Q = 1.4 \text{ nC}$
 $\sigma_x = 49 \mu\text{m}$,
 $\sigma_y = 8.5 \mu\text{m}$,
 $\sigma_z = 3.5 \mu\text{m}$
 Peak field
 $\sim 50 \text{ GV/m}$






N_2 gas sheet
 Density = 10^{17} m^{-3}
 Thickness = $150 \mu\text{m}$
 $Q_{\text{ions}} = 35 \text{ fC}$ ($> 10^5$ ions)

- Substantially larger yield
- Reduced densities
- flexibility
- Gas sheet tunability


WARP simulations: Impact and field ionization considered for different regimes

Ion Microscope - Design


- Transport and magnify generated ion beam to Multi-Channel Plate (MCP) detector
- Simulated fields and ion beam transport in CST
- Tested magnification imaging ability with test beam with random initial transverse velocity
- Design magnification: 8x






- Test data for surrogate model (WARP)
 - Variances in Q , σ_x , σ_y , σ_z , n_p
- Convolutional neural network (CNN)
- Robust performance of CNN
- Further enhancement when constraints enforced (e.g. from experimental measurements)



Sample beam





3D parameter space



Performance plots: surrogate model vs truth, and relative R^2 values



Gas Sheet Generation - Validation

- Bench tests
 - Fast recovery of vacuum after gas sheet operation (N_2)
 - Gas density profile agrees for various skimmer sizes




Vacuum Readings – 0.5 Hz

Gas Jet



1mm "Bulldozer" Collector



Gas Profile at Diagnostic

Commissioning Status

- UCLA Pegasus beamline
 - Validate results with well-diagnosed beams
 - Impact ionization dominates
 - Expected 9 fC charge (~55k ions)
 - Resolvable on MCP/phosphor
- Status
 - Installed on beamline, vacuum <E-9
 - Ion microscope commissioned on dedicated laser stand

Pegasus Beam

$$\sigma_x = \sigma_y = 40 \mu\text{m}$$

$$\sigma_z = 200 \mu\text{m}$$


$$Q = 0.1\text{nC}$$

Gas Sheet


$$n_{\text{N}_2} = 10^{14} \text{ cm}^{-3}$$

$$L_{\text{N}_2} = 150 \mu\text{m}$$



$$\theta = 45^\circ$$



WARP simulations for Pegasus parameters




GSM at Pegasus




Microscope Commissioning

- GSM is viable diagnostic
- Bench tested individual components
 - Skimmers / pumping
 - Ion microscope
 - Reconstruction methods
- First tests at Pegasus UCLA
- Next tests at SLAC FACET II
 - High intensity
- Opportunities for integration in beamline feedbacks with machine learning
 - Unique source of beam information



Initial Ion Capture on Phosphor Screen



Integrated Phosphor Intensity vs. Gas Density

N. Burger, G. Andonian, T. Hodgetts, D. Gavryushkin, J. Penney, A.L. Lamure, N. Norvell, M. Ruelas (RadiaBeam), A. Diaw, N. Cook, C. Hall (RadiaSoft), A. Marinelli, D. Cesar, M. Hogan, V. Yakimenko (SLAC), P. Denham, D. Garcia, P. Musumeci, A. Ody (UCLA), C. Welsch, M. Yadav (Univ. Liverpool)

UCLASLACradiasoft

UNIVERSITY OF
LIVERPOOL

This work supported by Grant No. **DE-SC0019717**